

Lab session 8:


Group-level SPM contrasts

Andrew Bauer

03/16/16

Session no.	Date (all Wednesday)	Topic/activity	Topic of quiz that day	Topic of lab write-up (assignment) due that day
1	13-Jan	Lab overview		
2	20-Jan	Brain anatomy		
3	27-Jan	Data preprocessing	Brain anatomy (no. 1)	
4	3-Feb	Set up GLM model	Functional brain anatomy (no. 2)	
5	10-Feb	Single-subject SPM contrasts	Data preprocessing and GLM model (no. 3)	Brain anatomy (no. 1)
6	17-Feb	Within-subject MVPA		Single-subject SPM contrasts (no. 2)
7	24-Feb	SIBR tour and review for mid-term exam		Within-subject MVPA (no. 3)
No lab	2-Mar	No lab (mid-term exam)		
No lab	9-Mar	No lab (spring break)		
8	16-Mar	Group-level SPM contrasts		
9	23-Mar	Between-subjects MVPA		Group-level SPM contrasts (no. 4)
10	30-Mar	Voxel-wise modeling		Between-subjects MVPA (no. 5)
11	6-Apr	Functional connectivity analysis (no assignment)		
12	13-Apr	Review for final exam		Voxel-wise modeling (no. 6)
No lab	20-Apr	No lab		
No lab	27-Apr	No lab (final exam)		

General sequence of data preprocessing and GLM analysis



Why analyze a group of subjects?

- The goal of a psychological research study is to evaluate a falsifiable hypothesis about a *population* (whose members vary in many respects)
 - A study is conducted on a *sample* of some population...
 - ... and a conclusion about that population is made by generalizing from the sample
 - The sample should be as representative of the population as possible

Getting a representative sample

- Can recruit subjects from different parts of the world using Amazon Mechanical Turk
 - Reference to late 18th century fake chess-playing machine (a person hid within the machine and did all the playing)

Mechanical Turk is a marketplace for work.

We give businesses and developers access to an on-demand, scalable workforce. Workers select from thousands of tasks and work whenever it's convenient.


240,837 HITs available. [View them now.](#)

Make Money by working on HITs

HITs - *Human Intelligence Tasks* - are individual tasks that you work on. [Find HITs now.](#)

As a Mechanical Turk Worker you:

- Can work from home
- Choose your own work hours
- Get paid for doing good work



Get Results from Mechanical Turk Workers

Ask workers to complete HITs - *Human Intelligence Tasks* - and get results using Mechanical Turk. [Register Now](#)


As a Mechanical Turk Requester you:

- Have access to a global, on-demand, 24 x 7 workforce
- Get thousands of HITs completed in minutes
- Pay only when you're satisfied with the results



How will we go about analyzing data on the group level?



(1) Determine how to *prepare* subjects' data that allows for group-level analysis



(2) Determine how to *treat* each subject's prepared data in the analysis

How to account for neuroanatomical differences between subjects? **Use spatial normalization**

- Morphs an individual's brain to a common template
 - We are using the Montreal Neurological Institute (MNI) template, which is the average of 152 brains
- Examine activation in the *same voxels* across subjects



Instead of spatial normalization, one could also identify the same brain regions in each subject using “functional localizers”

- Use a localizer scan to identify a brain region of interest in each subject (e.g. show pictures of faces to find each subject’s “face area”); the regions will be slightly different across subjects
- Then, extract activation data for different trials/conditions from this region in each subject
- *Advantage:* An individual’s data are left intact; no warping of the data into a template brain, which could introduce errors and noise
- *Disadvantages:* May be difficult to ensure that the localizer reveals valid brain areas (“is this *really* the face area for subject X?”); also, you don’t get nice voxel-level spatial resolution for comparing across subjects (as in spatial norm.)

Faces > Objects Passive Viewing





Image from Wojciulik, Kanwisher, & Driver (1998)

How will we go about analyzing data on the group level?

(1) Determine how to *prepare* subjects' data that allows for group-level analysis



(2) Determine how to *treat* each subject's prepared data in the analysis

Fixed-effects analysis (X now discouraged)

Subjects' fMRI
time courses

Averaged (combined)
time courses

Statistical map (from
averaged time courses)
(e.g. *Tools - Buildings*)




Image from Huettel, Song, & McCarthy: Functional Magnetic Resonance Imaging, 2nd ed.

- Averages subjects' data, giving equal weight to each subject; assumes that any effect is the same (“fixed”) for each subject, apart from some noise
- By averaging, this analysis does **NOT** account for the *variation* of the effect across subjects, which is a defining element of a *population*
 - Thus, this analysis is heavily influenced by unrepresentative outlier subjects

Random-effects analysis (✓ standard practice)

Subjects' fMRI
time courses

Subjects' statistical maps
(e.g. *Tools - Buildings*)

Combined (i.e. group-level)
statistical map
(e.g. *Tools - Buildings*)





Image from Huettel, Song, & McCarthy: Functional Magnetic Resonance Imaging, 2nd ed.


- Does **NOT** average subjects' data; treats pooled subjects as a *distribution*
- **DOES** account for the variation of the effect across subjects
 - This analysis recognizes outliers in the distribution and is less influenced by them

How we will analyze group-level data in the lab


(1)



Normalization




(2) *Random-effects analysis*
in each voxel (in MNI space)



What are group-level results generally like?

- Stronger activation, which often means:
 - Greater number of statistically significant voxels (using the same p -value threshold when comparing a single-subject to a group-level contrast)



- Also, group-level results contain *common* clusters of activation
 - **Disappearance of clusters that correspond to only one or a couple subjects**
 - The less common a cluster, the less likely it will be statistically significant in the group-level analysis

Multiple comparisons problem

(Applies to *both* single-subject and group-level analysis)

- Say we set $p = 0.05$ threshold for each voxel in a contrast
 - This means we are willing to tolerate a 5% chance that a voxel is statistically significant purely due to chance
 - If 200,000 voxels tested in an analysis, then there could be $200,000 * 0.05 = \mathbf{10,000}$ **falsely significant voxels**
- Ways to combat this problem:
 - Lower the p -value threshold to make a statistically significant voxel less likely due to chance
 - Set a cluster size threshold (e.g. 10, 20, 30 voxels)
 - A big cluster of significant voxels is less likely to be due to chance

Now follow the slides below to generate a group-level contrast as practice for the assignment...

Animals vs. other object categories
(group level, $N = 10$ subjects)

Start Matlab 2012b (on desktop, or type "matlab" in Start menu to find it)

NOTE: You MUST select Matlab 2012b, do NOT select 2014b

MATLAB 7.10.0 (R2010a)

File Edit Debug Parallel Desktop Window Help

Current Folder: C:\Users\andrewba\Documents\MATLAB

Shortcuts How to Add What's New

Command Window

New to MATLAB? Watch this [Video](#), see [Demos](#), or read [Getting Started](#).

MATLAB desktop keyboard shortcuts, such as Ctrl+S, are now customizable.
In addition, many keyboard shortcuts have changed for improved consistency across the desktop.

To customize keyboard shortcuts, use [Preferences](#). From there, you can also restore previous default settings by selecting "R2009a Windows Default Set" from the active settings drop-down list. For more information, see [Help](#).

[Click here](#) if you do not want to see this message again.

>>

Workspace

Name Value Min

Command History


```
%-- 2/21/13 7:08 PM --%
spm fmri
%-- 2/22/13 2:30 PM --%
1-tcdf(2.75,30)
1-tcdf(3,30)
1-tcdf(3,40)
1-tcdf(2.75,40)
%-- 2/25/13 10:30 AM --%
spm fmri
clc
pwd
ll
ls
clc
%-- 2/27/13 3:25 PM --%
```

Details

Start Ready


OVR

Select Set Path




1: Select Add with Subfolders

2: SINGLE-click the folder spm or spm8 under C:, click OK




1: Select Add with Subfolders again



1: Select Close

2: Select No if it asks to save the path file somewhere else (it may not ask)



1: Go to the Matlab Command Window and type:

`cd C:/Users/Public/lab_85_429_2016/session8`

...(OR navigate there using the browser)

2: Then type: `spm fmri`

The screenshot displays the MATLAB R2013a environment. The top toolbar includes options like 'New Script', 'Open', 'Find Files', 'Import Data', 'Save Workspace', 'Clear Workspace', 'Analyze Code', 'Run and Time', 'Clear Commands', 'Simulink Library', 'Layout', 'Set Path', 'Parallel', 'Help', 'Community', 'Request Support', and 'Add-Ons'. A tooltip 'Change the search path used by MATLAB to look for files' is visible near the 'Set Path' button. The 'Current Folder' pane on the left shows the path 'C:\Users\andrewba\Documents\MATLAB'. The 'Command Window' is active, showing a prompt 'New to MATLAB? Watch this Video, see Examples, or read Getting Started.' followed by the command `>> restoredefaultpath;matlabrc` and several subsequent '>>' prompts. The 'Workspace' pane on the right is empty. The 'Command History' pane at the bottom right shows a list of commands, including `delete modelSpec_wrongNoImage`, `spm fmri`, and `restoredefaultpath;matlabrc`, with timestamps for each execution. Red arrows point from the instructional text to the Command Window and the Command History pane.

NOTE: When typing commands, press the Tab key for auto-completion (may sometimes only partially complete what you're typing)

1: Select
Specify 2nd_
level

2: Select Directory

3: Select animal_0013,
then press Done

4: Next
slide...

The image is a screenshot of the SPM8 software interface, showing several windows and their contents. The windows include:

- SPM8 (andrewba): Menu**: A window with various processing options like Realign, Slice timing, Smooth, Coregister, Normalise, Segment, Specify 1st-level, Review, Specify 2nd-level, Estimate, Results, and Dynamic Causal Modelling.
- Batch Editor**: A window with a menu bar (File, Edit, View, SPM, Basic) and a Module List. The current module is 'Factorial design specification'. The 'Scans' directory is selected in the list.
- File Explorer**: A window showing the directory structure of 'E:\lab_85_429_2014\'. It shows subdirectories like '04383B' and '04408B'. A search filter '^con_0013.*' is applied, and 11 matching files are listed.
- Batch Editor - Factorial design specification**: A window showing the 'Current Module: Factorial design specification' and a list of options: Directory, Design, One-sample t-test, Scans, Covariates, Masking, Threshold masking, None, Implicit Mask, Explicit Mask, Global calculation, Omit. The 'Current Item: Scans' is highlighted.

Annotations and steps are provided in red boxes with arrows pointing to the relevant parts of the interface:

- Step 5: Select Scans**: Points to the 'Scans' directory in the Batch Editor's Module List.
- Step 6: Select the ".." (two periods) to go up a directory**: Points to the '..' directory in the File Explorer.
- NOTE: Leave as 1**: Points to the '1' in the search filter '^con_0013.*'.
- Step 7: Type: ^con_0013.* then press Rec, then press Done (all 10 subject contrast images will be selected for you)**: Points to the search filter '^con_0013.*' and the 'Rec' and 'Done' buttons in the File Explorer.
- Step 8: Click Play, then exit Batch Editor, don't save (here or ever). NOTE: if it asks you to overwrite, click "continue" to overwrite**: Points to the 'Play' button in the Batch Editor's menu bar.
- NOTE: the contrast number 13 matches the category number for Animals, which is 13**: Points to the '13' in the search filter '^con_0013.*'.


1: After the group-level contrast is done being specified (previous slide, it takes only a second or two), select Estimate

5: Click Play

2: Select SPM.mat

3: Select animal_0013 folder if you're not already there (NOTE: Your screen should show session8, not session9)

4: Then select SPM.mat, then press Done



1: After the group-level contrast is finished being estimated (previous slide), select Results


2: Select animal_0013 folder if you're not already in that directory

3: Select SPM.mat, then press Done

4: Next slide...




5: Select Define new contrast



6: Type the name for the contrast (after the "define contrast" window appears): animals_vs_others

7: For the contrast value, type: 1 then press OK






NOTE: Make sure that the t-contrast button is selected

8: Next slide...

The image shows a screenshot of the SPM8 (eprime) software interface. The main window is titled "SPM8 (eprime): Graphics" and contains a menu bar with File, Edit, View, Insert, Tools, Desktop, Window, SPM Figure, and Help. Below the menu bar is a large empty area. In the foreground, there is a smaller window titled "SPM contrast manager" with a subtitle "Select contrasts...". This window has a tabbed interface with "t-contrasts", "F-contrasts", and "all" tabs. The "t-contrasts" tab is selected, showing a list of contrasts. The first entry, "001 (T) : animal_vs_other", is highlighted. A red arrow points from this entry to a red box at the bottom of the slide. To the right of the list is a "contrast(s)" field, a "Design matrix" field, and a "parameter estimability" field. At the bottom of the dialog are buttons for "Define new contrast...", "Reset", and "Done". A status bar at the bottom of the dialog says "Selected 1 contrast, press 'Done' when finished." In the background, other SPM8 windows are visible, including "SPM8 (eprime): Menu" and "SPM8 (eprime): Stats: Results".

1: Click the animals_vs_others entry so that it's highlighted; then press Done

2: Next slide...



3: Apply masking:
none, then enter

4: Title for comparison: keep at
default "animals_vs_others"

5: Threshold: First choose "none" for " p value adjustment to control", then type: 1

(we will view results in xjview and set the p threshold there)

6: Extent threshold: 0
(we will set it in xjview)

1: After defining the contrast parameters in the previous slide, spm will display a window of results. However, we will ignore this and instead view results using xjview, as we did with single-subject spm contrasts. But don't close spm. Go back to the Matlab command window (it's already open) and cd to the directory with the results. That is, type:

```
cd C:/Users/Public/lab_85_429_2016/session8/animal_0013
```

(Unless you're already in this folder, in which case, you're fine)

The screenshot shows the MATLAB environment. The Command Window displays the following commands and output:

```
>> pwd  
ans =  
E:\lab_85_429_2014\session8\animal_0013  
  
>> ls  
.  
..  
RPV.hdr  
RPV.img  
SPM.mat  
spmT_0001.hdr  
spmT_0001.img  
RPV.img  
ResMS.hdr  
ResMS.img  
beta_0001.hdr  
beta_0001.img  
mask.hdr  
con_0001.hdr  
con_0001.img  
mask.hdr  
spmT_0001.hdr  
spmT_0001.img  
mask.hdr  
spmT_0001.hdr  
spmT_0001.img  
  
>> xjview  
fx >>
```


Annotations in the image include:

- A red circle around the directory path in the Command Window: `E:\lab_85_429_2014\session8\animal_0013`.
- A red circle around the `pwd` command.
- A red circle around the `ls` command.
- A red circle around the `xjview` command.

2: To see what directory you're in, look at this red-circled area, OR type: `pwd`
(NOTE: your correct directory will be slightly different from the example directory shown here in this slide)

3: To see what files are in this directory, type: `ls`

4: If the directory is correct and you see files such as "spmT_0001.img", then next type: `xjview`



1: After xjview starts up, to start viewing results, go to File->Open image

NOTE: Leave as 1

2: Select spmT_0001.img and press Done
(This is the contrast image that we created using spm; it is a image of *t* scores after a voxel-wise *t*-test (over all subjects) between Animals and the average of the other categories)


NOTE: For any of these group-level category contrasts (animal, tool, etc.), you will calculate *one* contrast *only* and the result will *always* be saved as spmT_0001.img into the corresponding category folder. If you made a mistake in the first contrast, then you will have to define it again and use spmT_0002.img *instead*

3: Next slide...

xjView: E:\lab_85_429_2014\session8\animal_0013\spmT_0001.img.1

File Edit View Insert Tools Desktop Window Help

1: Check Render View and select "old" (instead of "new") right next to the check




// Left Cerebrum // Occipital Lobe // Sub-Gyral // White Matter // undefined // **Fusiform_L (aal)**

XHairs Off
avg152T1
avg152T2
avg305T1
ch2
ch2bet
aal
brodmann
other...
colorbar max
auto

E:\lab_85_429_2014\session8\animal_0013\spmT_0001.img.1
This is a T test image.
mat =
-2 0 0 80
0 2 0 -114
0 0 2 -52
0 0 0 1
dimension =
79 95 69

search Sub-Gyral in xBrain
overlay Amygdala
report volume Small volume common region slice view
display intensity ☐ All ☒ Only + ☐ Only - ☒ Render View **old**
cluster size 10 Pick Cluster/... Select Cluster Clear Selection
pValue: 0.005 FDR p= T= 3.2498 df= 9
x = -39.43 y = -56.10 z = -14.97





NOTE: The window won't maximize on our computers

NOTE: Wherever you click in the brain (in one of the three slice views), the MNI coordinates will adjust themselves accordingly. You can also manually enter coordinates into these fields and your cursor will go to that point (although the coordinates may change slightly, but that's okay – don't worry about it)

NOTE: The AAL brain region name will appear in this circled area whenever your cursor is within the grey matter of the brain. Sometimes the AAL name won't appear; it means that that point isn't defined there, and your cursor is probably in white matter or on/outside the skull

2: Next slide...





1: To show a single cluster and display its anatomical information, first click within *any* cluster on the screen

2: Click Select Cluster. It should say that there is one cluster selected in the white space to the left. (If you must de-select a cluster ever, click Clear Selection to the right)

3: Click Pick Cluster/Info. Now you should see # voxels per brain area in the white space to the left


4: Pay attention *only* to the AAL brain areas and the TOTAL # VOXELS value

Important notes for assignment

- *Always...*
 - Load the file `spmT_0001.img` in xjview for a given category group contrast, found in that category folder
 - Unless you made a mistake and had to create the contrast a second time, in which case use `spmT_0002.img`, etc.
 - Display *positive* activation (toggle the “Only +” button)
 - Use a cluster size threshold of 10 voxels
 - Unless stated otherwise
 - Use a p -value threshold of 0.005
 - Unless stated otherwise

- This is all the guidance that you should need to do this session's assignment
- If you haven't already, you are free to download the xjview manual off Blackboard for further help
- If something goes wrong with spm *OR* xjview, just exit the program and restart it from Matlab as usual. You don't have to exit out of Matlab to restart spm or xjview
- If you must restart Matlab for whatever reason, after you start Matlab up *be sure to set the paths again and cd to the directory of whatever category contrast you are working on*
- See the slides below for help with orientation terms, planes, brain anatomy, and general functions of brain regions


Orientation terms and planes




We are using the **Automated Anatomical Labeling (AAL)** atlas in this lab (very similar to the atlases of your two quizzes)



Lateral gyri and sulci



Medial gyri (some redundancy w/previous slide)



General functional neuroanatomy

