

# Lab session 4:

## Set up General Linear Model (GLM)

Andrew Bauer

02/03/16


Session no.	Date (all Wednesday)	Topic/activity	Topic of quiz that day	Topic of lab write-up (assignment) due that day
1	13-Jan	Lab overview		
2	20-Jan	Brain anatomy		
3	27-Jan	Data preprocessing	Brain anatomy (no. 1)	
4	3-Feb	Set up GLM model	Functional brain anatomy (no. 2)	
5	10-Feb	Single-subject SPM contrasts	Data preprocessing and GLM model (no. 3)	Brain anatomy (no. 1)
6	17-Feb	Within-subject MVPA		Single-subject SPM contrasts (no. 2)
7	24-Feb	SIBR tour and review for mid-term exam	Study slides 6, 19, 20 below	Within-subject MVPA (no. 3)
No lab	2-Mar	No lab (mid-term exam)		
No lab	9-Mar	No lab (spring break)		
8	16-Mar	Group-level SPM contrasts		
9	23-Mar	Between-subjects MVPA		Group-level SPM contrasts (no. 4)
10	30-Mar	Voxel-wise modeling		Between-subjects MVPA (no. 5)
11	6-Apr	Functional connectivity analysis (no assignment)		
12	13-Apr	Review for final exam		Voxel-wise modeling (no. 6)
No lab	20-Apr	No lab		
No lab	27-Apr	No lab (final exam)		

## Data preprocessing pipeline




Generalized  
Linear  
Model (GLM)

Setting up analysis



Data analysis & results

# Our analysis: Discover the brain areas that are associated with thinking about different object concept *categories*




Instructions to subjects:  
“Think of the main properties of each named object”

**Table 1.** 60 stimulus words grouped into 12 semantic categories.

Category	Exemplar 1	Exemplar 2	Exemplar 3	Exemplar 4	Exemplar 5
body parts	leg	arm	eye	foot	hand
furniture	chair	table	bed	desk	dresser
vehicles	car	airplane	train	truck	bicycle
animals	horse	dog	bear	cow	cat
kitchen utensils	glass	knife	bottle	cup	spoon
tools	chisel	hammer	screwdriver	pliers	saw
buildings	apartment	barn	house	church	igloo
building parts	window	door	chimney	closet	arch
clothing	coat	dress	shirt	skirt	pants
insects	fly	ant	bee	butterfly	beetle
vegetables	lettuce	tomato	carrot	corn	celery
man-made objects	refrigerator	key	telephone	watch	bell


# The *H*emodynamic *R*esponse *F*unction of the BOLD MRI signal

- MRI signal in the brain increases by a few % points after there are increases in neural activity



# Fitting the HRF to actual MRI signal (Quiz no. 3 slides 6, 19, 20)

- If a voxel is active, we expect its MRI signal to take the *same* HRF shape, regardless of the brain region it's in or which subject's data it is
- “Brain activation level” of condition  $x$  in voxel  $i$  = how well the HRF *fits* (and the magnitude/height of the fitted waveforms) to the actual MRI signal during that condition's trials



# What is “brain activation level” in fMRI?

- Voxel-wise GLM regression model:  $y = \beta_1 x_1 + \beta_2 x_2 + \beta_3 x_3 + \dots + \beta_n x_n + \varepsilon$

$y$ : MRI signal (time series) at a voxel – **different** across voxels

$x$ : independent variables/conditions (regressors) – **same** across voxels

$\beta$ : regression coefficients (HRF fits/magnitudes) – **different** across voxels

$\varepsilon$ : residuals (data unaccounted for by model) – **different** across voxels

We are measuring this

- So the  $\beta$  coefficient is the ***brain activation level*** associated with a particular condition  $x$

- In our case, the conditions are the different object concept *categories*
- The model estimates the  $\beta$  value for each condition by how well the HRF fits (and its magnitude) the MRI signal over all a given condition's trials

Start Matlab 2012b (on desktop, or type "matlab" in Start menu to find it)

**NOTE: You MUST select Matlab 2012b, do NOT select 2014b**

MATLAB 7.10.0 (R2010a)

File Edit Debug Parallel Desktop Window Help

Current Folder: C:\Users\andrewba\Documents\MATLAB

Shortcuts How to Add What's New

Command Window

New to MATLAB? Watch this [Video](#), see [Demos](#), or read [Getting Started](#).

MATLAB desktop keyboard shortcuts, such as Ctrl+S, are now customizable.  
In addition, many keyboard shortcuts have changed for improved consistency across the desktop.

To customize keyboard shortcuts, use [Preferences](#). From there, you can also restore previous default settings by selecting "R2009a Windows Default Set" from the active settings drop-down list. For more information, see [Help](#).

[Click here](#) if you do not want to see this message again.

>>

Workspace

Stack: Select data to plot

Name Value Min

Command History

```
%-- 2/21/13 7:08 PM --%
spm fmri
%-- 2/22/13 2:30 PM --%
1-tcdf(2.75,30)
1-tcdf(3,30)
1-tcdf(3,40)
1-tcdf(2.75,40)
%-- 2/25/13 10:30 AM --%
spm fmri
clc
pwd
ll
ls
clc
%-- 2/27/13 3:25 PM --%
```


Details

Start Ready

OVR




Select Set Path




1: Select Add with Subfolders

2: SINGLE-click the folder spm or spm8 under C:, click OK




1: Select Add with Subfolders again



1: Select Close


2: Select No if it asks to save the path file somewhere else (it may not ask)



```
cd C:/Users/Public/lab_85_429_2016/session4
```

2: Then type: `spm_fmri`

...(OR navigate there using the browser)



1: Select  
Specify 1st-  
level

2: Select Directory


4: For Units for  
design, select Seconds

5: For interscan  
interval (TR),  
type: 1

3: Select 04647B, then  
press Done

6: Next  
slide...





The screenshot shows the SPM8 software interface with three windows open:


- SPM8 (andrewba): Menu**: Contains options like Realign, Slice timing, Smooth, Coregister, Normalise, Segment, Specify 1st-level, Review, Specify 2nd-level, Estimate, Results, and Dynamic Causal Modelling.
- Batch Editor**: Shows a Module List with 'fMRI model specification' selected. A red arrow points to the 'Play' button (green triangle) in the top toolbar.
- Multiple conditions**: A file explorer window showing the directory 'C:\Users\andrewba\Desktop\85\_429\_2014\session4\'. It lists files like '04647B' and 'conditionInfo.mat'. A red arrow points to the 'Done' button at the bottom.

Three red text boxes with arrows provide instructions:

- Box 13** (points to the Play button): "13: Click Play, then exit Batch Editor, don't save (here or ever), then you're done... a design matrix window should appear after some time"
- Box 11** (points to 'Multiple conditions' in the Module List): "11: Select Multiple conditions"
- Box 12** (points to the 'conditionInfo.mat' file): "12: Select conditionInfo.mat, then press Done"



## Statistical analysis: Design



·Sn(1) fixation*bf(1)	Baseline
·Sn(1) animal*bf(1)	Animal
·Sn(1) bodypart*bf(1)	Bodypart
·Sn(1) building*bf(1)	Building
·Sn(1) buildingpart*bf(1)	Buildingpart
·Sn(1) clothing*bf(1)	Clothing
·Sn(1) furniture*bf(1)	Furniture
·Sn(1) insect*bf(1)	Insect
·Sn(1) kitchen*bf(1)	Kitchen
·Sn(1) manmade*bf(1)	Manmade
·Sn(1) tool*bf(1)	Tool
·Sn(1) vegetable*bf(1)	Vegetable
·Sn(1) vehicle*bf(1)	Vehicle
·Sn(1) constant	Constant

# Programming statistical $\beta$ contrasts to compare...

- Activation of one condition vs. Baseline (“Resting”)
- Activation between two different conditions
  - Baseline cancels out (see below)

Contrast	Baseline multiplier	Tools multiplier	Buildings multiplier	Statistical result ( <i>t</i> -stat and <i>p</i> -value) in <u>each</u> voxel
$\beta_{\text{Tools}} - \beta_{\text{Baseline}}$	-1	1	0	<b>Tools activation that is &gt; “Resting”</b>
$\beta_{\text{Buildings}} - \beta_{\text{Baseline}}$	-1	0	1	<b>Buildings activation that is &gt; “Resting”</b>
$\beta_{\text{Tools}} - \beta_{\text{Buildings}}$ i.e. (Tools - Baseline) - (Buildings - Baseline)	0	1	-1	(Tools > “Resting”) that is > (Buildings > “Resting”)
$\beta_{\text{Buildings}} - \beta_{\text{Tools}}$ i.e. (Buildings - Baseline) - (Tools - Baseline)	0	-1	1	(Buildings > “Resting”) that is > (Tools > “Resting”)

# Data preprocessing sequence (Quiz no. 3

slides 6, 19, 20)

## 1. Motion correction

- Ensures that brain activity of voxel  $i$  corresponds to the same volume of brain tissue throughout whole experiment

## 2. Slice time correction

- Need to correct for fact that brain slices are not collected simultaneously each image acquisition (TR)

## 3. Spatial normalization

- Morphs each subject's brain to a common/template brain to examine activation in same locations across subjects

## 4. Smoothing

- Reduces amount of noise per voxel by averaging over nearby voxels

# Final notes about data preprocessing

(Quiz no. 3: Slides 6, 19, 20)

- In brain imaging and other fields, data processing can result in ***loss*** or ***artificial creation*** of data
  - Smoothing (you control how much)
    - ***Lose*** small fluctuations in data that might be real
    - If one voxel's activation level is high due to error/random chance, then averaging it in with other voxels will spread this error to other voxels, ***creating an artificial*** cluster of activation
  - Spatial normalization
    - ***Lose*** data due to the normalization algorithm tampering with the original data
    - ***Lose*** individual quirks of a person's data

# References

Just, M. A., Cherkassky, V. L., Aryal, S., & Mitchell, T. M. (2010). A neurosemantic theory of concrete noun representation based on the underlying brain codes. *PLoS One*, 5(1), e8622.  
doi:10.1371/journal.pone.0008622